MEMORANDUM

date: 2020 02 14

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Please note that there is a new test that is being introduced for TB PCR, as of February 24, 2020. This new test is called the GeneXpert and it replaces the Roche TB Amplicor test.

Although the machine is different, the test is offered for the same two indications as before: 1) sputum samples that are positive for acid-fast bacilli on microscopic examination, and; 2) sputum samples from a patient who is considered by the treating physician to be high-risk for active TB.

There are three differences in performance that clients might note. First, the test has a higher analytic sensitivity than the old Roche-based test. Therefore, someone who was previously PCR- might be PCR+ with the new test. This does not mean the patient is sicker or more contagious. It just means that the new test can detect a smaller number of bacteria in a sample than the previous test. Second, the new test can be done on a single sample, instead of waiting for a week's worth of samples to be processed in a batch. You might note that instead of waiting 3-5 days for a PCR, you can expect to get results within one working day of the microscopy result. Third, the new test both looks for *Mycobacterium tuberculosis* DNA and reports if there is a mutation in a gene called RNA polymerase B (rpoB) that is associated with rifampin resistance.

As a result, the report will be modified to have two parts: a) *M. tuberculosis* detected (or not detected) and b) Rifampin resistance detected (or not detected).

The new test has been validated on sputum samples that are positive for acid-fast bacilli that subsequently grew *M. tuberculosis* (positive controls) or subsequently grew non-tuberculous mycobacteria (negative controls). It has also been validated to detect the rifampin resistance causing mutation. The test has not been validated on other samples and is subject to the same limitations observed with other PCR assays in that sensitivity is compromised in paucibacillary infections and extrapulmonary samples.

If you have any questions or concerns, please do not hesitate to contact the Clinical Laboratories Client Services department at 514-934-1934, ext 35687.