## Génotypage *DPYD* (4 variations récurrentes) – INSTRUCTIONS

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## *DPYD* Genotyping (4 recurrent variants) – INSTRUCTIONS

The MSSS *Direction de la biovigilance et de la biologie médicale* has mandated the MUHC (Glen site) Core Molecular Diagnostic Laboratory (CMDL) to offer *DPYD* genotyping (4 recurrent variants) for detection of severe fluoropyrimidine toxicity due to dihydropyrimidine dehydrogenase (DPD) deficiency.

**Please read carefully the specifications and limitations of this test, provided below.**

**Please follow the provided instructions, not to delay results and treatment.**

For more information, please contact the persons responsible for this test at the CMDL:

- **Andrea Gomez**, Clinical biology specialist  
[andrea.gomez@muhc.mcgill.ca](mailto:andrea.gomez@muhc.mcgill.ca), 514-934-1934 x38776
- **William Foulkes**, Medical geneticist  
[william.foulkes@mcgill.ca](mailto:william.foulkes@mcgill.ca), 514-934-1934 x44067

Thank you for your collaboration. **Please forward this information to all relevant stakeholders.**

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## DPYD Genotyping (4 recurrent variants) – INSTRUCTIONS

Please read carefully the specifications and limitations of this test, provided below.  
Please follow the provided instructions, not to delay results and treatment.

|                          |   |
|--------------------------|---|
| <b>Description:</b>      | DPYD genotyping (4 recurrent variants)      |
| <b>Sector:</b>           | Core Molecular Diagnostic Laboratory (CMDL) |
| <b>Requisition form:</b> | Attached below                              |
| <b>MSSS code:</b>        | <a href="#">65036</a>                       |
| <b>LIS code:</b>         | MolGenDPYD                                  |
| <b>Turnaround time:</b>  | 5 working days, from sample reception       |

### TEST BACKGROUND

This test consists in genotyping 4 recurrent *DPYD* variants associated with an increased risk of severe fluoropyrimidine toxicity, including:

- 2 variants (c.1905+1G>A/DPYD\*2A and c.1679T>G) with no DPD function (activity value = 0).
- 2 variants (c.2846A>T and c.1236G>A, HapB3 proxy) with decreased DPD function (activity value = 0.5).

The goal of this test is to adjust fluoropyrimidine treatment in individuals with one or more variants, depending on the *DPYD* genotype. The INESSS provides a [clinical tool](#) for information purposes

### TEST METHOD

The 4 variants are genotyped by a TaqMan assay in duplicate (2 independent DNA extractions, hence the importance of providing **2 tubes of blood**) from 4 independent PCR reactions (one for each variant). Results (PDF report) are faxed to the **fax number provided by the referring physician on the requisition form** (see form below).

### LIMITATIONS

This test is known to have several important limitations, notably:

**A limited sensitivity**, regardless of the ethnic origin of the person tested: Over 100 *DPYD* variants known to alter DPD activity have been described to date and other functional variants remain to be discovered. Furthermore, fluoropyrimidine metabolism, efficacy and risk of toxicity may be influenced by genetic and non-genetic factors that are not assessed by this test.

**A significant ethnic bias**: This test detects the 4 most frequent *DPYD* drug response variants in individuals of European descent and does not detect some of the most frequent variants in other ethnic groups, notably in Africans and East Asians. This test has limited to no clinical validity in individuals of non-European descent, depending on the ethnic origin.

**A risk of rare false-positive results**: The only patients with complete DPD deficiency potentially detected by this genotyping method are individuals homozygous for one of the 2 variants with no DPD function (c.1905+1G>A/DPYD\*2A and c.1679T>G). For individuals with 2 heterozygous variants, genotyping alone cannot determine if these variants are on different alleles (in trans,

complete DPD deficiency with no activity) or on the same allele (in cis, partial DPD deficiency, with 50% residual activity).

### INCIDENTAL FINDINGS

Individuals with 1 heterozygous null *DPYD* allele are unaffected heterozygous carriers of a rare inborn error of metabolism (dihydropyrimidine dehydrogenase deficiency, autosomal recessive, [MIM 274270](#)). Carrier testing could be recommended for at-risk relatives and reproductive partners of known carriers, along with **genetic counselling**.

### SAMPLES & SHIPPING

**Sample:** Blood draw

**Regular container:** **2 tubes** (EDTA, purple top)

**Required volume:** 2 x 4 mL

**Fasting / special diet:** No

**Shipping delay:** Immediately

**Storage instructions:**

- Please send samples to the laboratory at room temperature as soon as possible.
- Do not centrifuge blood samples.
- Do not store at extreme temperatures (do not freeze or store at temperature > 30°C).

**Shipping address:**

Core Molecular Diagnostic Laboratory (CMDL), Glen site

OPTILAB-MUHC Genetics

1001 Decarie boul., E05.5051

Montreal, QC, H4A 3J1 Canada

**Rejection criteria:**

- No/incomplete requisition form (**all fields marked with an asterisk are mandatory**).
- Specimen tubes not labelled with at least **two identifiers**.