Technology Assessment Unit (TAU) of
the McGill University Health Centre
(MUHC)

18 November 2024

Use of N95 respirators past their shelf life

Health Technology Assessment
Brief Report
Report no. 94

Report prepared for the Technology Assessment Unit (TAU) of the McGill University Health Centre (MUHC)

by

Nisha Almeida, PhD

Mission Statement

The MUHC Health Technology Assessment Unit (TAU) advises hospital administrators and clinical teams in difficult resource allocation decisions. Using an approach based on independent, critical evaluations of the available scientific evidence and a transparent, fair decision-making process, novel and existing medical equipment, drugs and procedures used by healthcare professionals are prioritized on a continuous basis ensuring the best care for life with the best use of resources.

Brief Reports

Brief reports are prepared in response to urgent requests for information or to update previous reports with new evidence; in such cases an in-depth evaluation is not possible or unnecessary. They are reviewed by the Manager of TAU and the Chair of the Policy Committee, but are not submitted to or approved by the Policy Committee, and usually contain no recommendations.

Declaration of Conflicts of Interest

Members of TAU's research staff and policy committee declare no conflicts of interest.

Suggested Citation

Almeida N. Use of N95 respirators past their shelf life. Montreal (Canada): Technology Assessment Unit (TAU) of the McGill University Health Centre (MUHC); Nov 18, 2023. Report no. 94. 27 pages

ACKNOWLEDGEMENTS

The expert assistance of the following individuals is gratefully acknowledged:

- Kim Tanguay, N. M.Ed., Nursing Practice Consultant, Nursing Directorate, McGill University Health Centre

REPORT REQUESTOR

This report was requested by Kim Tanguay, Nursing Practice Consultant, Nursing Directorate on November 10, 2023.

TYPES OF RECOMMENDATIONS ISSUED BY THE TAU COMMITTEE

Type of recommendation	Explanation
Approved	<ul style="list-style-type: none"> Evidence for relevant decision criteria, including efficacy, safety, and cost, as well as context-specific factors such as feasibility, is sufficiently strong to justify a recommendation that the technology be accepted, used and funded through the institutional operating budget
Approved for evaluation	<ul style="list-style-type: none"> There is a reasonable <i>probability</i> that relevant decision criteria, including efficacy, safety, and cost, as well as context-specific factors such as feasibility, are favorable but the evidence is not yet sufficiently strong to support a recommendation for permanent and routine approval. The evidence is sufficiently strong to recommend a <i>temporary</i> approval in a restricted population for the purposes of evaluation, funded through the institutional operating budget.
Not approved	<ul style="list-style-type: none"> There is insufficient evidence for the relevant decision criteria, including efficacy, safety, and cost; The costs of any use of the technology (e.g. for research purposes) should not normally be covered by the institutional budget.

DISCLAIMER

The Technology Assessment Unit ("TAU") of the McGill University Health Centre ("MUHC") was created in order to prepare accurate and trustworthy evidence to inform decision-making and when necessary to make policy recommendations based on this evidence. The objective of the TAU is to advise the hospitals in difficult resource allocation decisions, using an approach based on sound, scientific technology assessments and a transparent, fair decision-making process. Consistent with its role within a university health centre, it publishes its research when appropriate, and contributes to the training of personnel in the field of health technology assessment.

The information contained in this report may include, but is not limited to, existing public literature, studies, materials, and other information and documentation available to the MUHC at the time it was prepared, and it was guided by expert input and advice throughout its preparation. The information in this report should not be used as a substitute for professional medical advice, assessment and evaluation. While MUHC has taken care in the preparation of this report to ensure that its contents are accurate, complete, and up-to-date, MUHC does not make any guarantee to that effect. MUHC is not responsible for any liability whatsoever, errors or omissions or injury, loss, or damage arising from or as a result of the use (or misuse) of any information contained in or implied by the information in this report.

We encourage our readers to seek and consult with qualified health care professionals for answers to their personal medical questions. Usage of any links or websites in the report does not imply recommendations or endorsements of products or services.

TABLE OF CONTENTS

Acknowledgements	i
Report Requestor	i
Types of Recommendations Issued by the TAU committee.....	ii
Disclaimer	ii
Table of Contents	iii
List of Figures.....	iv
Summary	v
Sommaire	ix
List of Abbreviations.....	xiii
1. Background.....	1
1.1 Reason for Brief Report	1
2. Policy and Evaluation Questions	2
2.1 Policy Question	2
2.2 Evaluation Questions (Objective of this report)	2
3. Methods	2
3.1 Literature search.....	2
4. Evidence from Scientific Literature	2
4.1 Results of the literature search	2
4.2 Factors that determine N95 respirator performance.....	3
4.3 Summary of NIOSH findings.....	4
5. Solutions for stockpiled respirators.....	5
5.1 NIOSH recommendations for stockpiled respirators	5
5.2 A framework for evaluating and testing stockpiled FFRs	5
6. Guidance for expired N95 respirators	6
6.1 Health Canada guidance ⁷	6
6.2 WHO guidance	7
7. Conclusions.....	7
8. Recommendations.....	8
References.....	10
Appendix	11

LIST OF FIGURES

Figure 1: Hypothetical stratified random sampling design for N95 FFRs 11

Figure 2. Example Respirator Visual Inspection Checklist 12

SUMMARY

BACKGROUND

The COVID-19 pandemic necessitated the stockpiling of personal protective equipment (PPE) as a way to ensure the availability of necessary quantities during the public health emergency. As a consequence, healthcare establishments now have large stockpiled reserves of PPE, including N95 respirators, that have since exceeded their manufacture-designated shelf life.

The MUHC currently has a stockpile of 283,100 N95 Chompak F550 respirators that have a manufacture-designated shelf life of 3 years, and that have recently expired or are due to expire in 2024. These masks were purchased by the MUHC at a cost of \$228,654.

Given the need to balance the safety and effectiveness of expired N95 masks against the considerable procurement and storage costs of PPE, it is necessary to evaluate whether N95 respirators past their shelf life continue to be effective in protecting healthcare personnel.

This evaluation was therefore requested by Kim Tanguay, Nursing Practice Consultant, Nursing Directorate of the McGill University Health Centre (MUHC) on November 10, 2023.

OBJECTIVES

The objectives of this report are to:

1. Identify the evidence on the effectiveness of expired N95 respirators in terms of aerosol filtration efficiencies;
2. Identify guidance for the safe use of N95 respirators past their expiry date.

FINDINGS

Filtration performance of expired N95 respirators

- The best available evidence on the viability of expired N95 respirators comes from a large study (n=3,971 stockpiled N95 respirators manufactured between 2003 and 2013) conducted by The National Institute for Occupational Safety and Health (NIOSH) in the US. They reported that:
 - 98% (3,895 of 3,971) of the evaluated respirators **maintained their inhalation and exhalation resistance and filtration performance**, with the upper limit of the 99% confidence intervals well within the NIOSH minimum performance requirements. Filtration performance was based on Standard Test Procedures

that designate pass/fail criteria for Particulate Filter Efficiency for N95 to be $\geq 95\%$, i.e. N95 FFRs must have a particle penetration of less than 5.0%.

- Storage conditions such as exposure to dust, moisture, sunlight and damage to pallets and cases are important factors that may impact respirator viability and should be carefully documented.
- While these results are suggestive that the majority of evaluated N95 respirators stockpiled for long periods of time maintain their protective properties, NIOSH cautions that their results may not be applicable to other respirator models not included in this study or to stockpile facilities with various storage conditions. **The Champak F550 N95 respirators were not included in this study.**

Guidance

- NIOSH acknowledges that managing stockpiled and expired respirators will be an ongoing consideration and that stockpile personnel require uniform guidance. They recommend:
 - The development of **an economically sustainable stockpile quality assurance process**, such that stockpiles from around the nation test and evaluate their respirator inventories using NIOSH's minimum performance requirements, and share their results with product manufacturers to optimize shelf life designations.
 - A national stockpile maintenance approach to purchase, store, and deploy respirators (e.g., rotate inventory) to ensure respirator availability and product usage within these optimized shelf lives.
- NIOSH developed a 4-step framework to guide stockpile managers to assess the performance of filtering facepiece respirators (FFR):
 1. Define the population of FFRs to be sampled;
 2. Use random sampling strategies to obtain a representative sample;
 3. Inspect and test the sampled units using standard testing procedures in accordance with government or Occupational Health and Safety norms and regulations;
 4. Compute results and extrapolate to the entire population to determine the probability of observing a similar failure rate in the whole population of FFRs.

CONCLUSIONS

- The best available evidence on the performance of expired or stockpiled N95 respirators is derived from a NIOSH evaluation of 3,971 stockpiled N95 respirators manufactured between 2003 and 2013, which found that 98% of the evaluated

respirators maintained their inhalation and exhalation resistance and filtration performance after their expiry date.

- Nonetheless, NIOSH cautions these results may not be applicable to other respirator models not included in this study or stockpile facilities with various storage conditions. The Champak F550 N95 respirators were not included in this study.
- Most regulatory bodies including Health Canada, the CDC and the WHO agree that respirators beyond their labelled shelf life, while no longer considered approved, may retain adequate filter performance if stored properly. However, because loss of elasticity in the bands or straps can impact a proper seal, they recommend that respirators be closely inspected for signs of damage and fit-tested before use.
- To guide stockpile personnel, NIOSH recommends the creation of a quality assurance process for stockpiled respirators to determine optimal shelf life designations.
- NIOSH has also developed a 4-step framework to guide stockpile managers in evaluating the performance of stockpiled respirators, such that a representative sample of respirators can undergo quality control to determine the probability of observing a similar failure rate in the whole population of respirators.

KEY FINDINGS AND RECOMMENDATIONS

- There is reasonable evidence to indicate that N95 respirators past their designated shelf life **retain adequate filter performance** if stored properly. NIOSH reported that 98% of almost 4000 evaluated N95 respirators past their expiry date maintained their protective properties. However, these results may not be generalizable to respirators not included in their analysis, including Champak F550 N95 respirators.
- **Recommendations for the MUHC:**
 - Stockpile managers should perform **quality assurance testing** on representative samples of Champak F550 N95 respirators from different lots and storage conditions, as described by NIOSH, so the viability of the complete stockpile can be ascertained. Representative sample testing allows extrapolation of findings to the population of respirators from which the sample was derived, allowing for a cost-effective solution;
 - The following factors should be evaluated:
 - Visual signs of damage to respirators and fit-testing before use
 - Filtration efficiency
 - Storage conditions
- **Recommendations at the national/province level (NIOSH recommendation):**
 - The creation of **an economically sustainable stockpile quality assurance process**, to ensure that stockpiles of respirator inventories across the province are evaluated using NIOSH's minimum performance requirements. Results should be shared with product manufacturers to optimize shelf life designations.
 - A stockpile maintenance approach at the national level to purchase, store, and deploy respirators (e.g., by rotating inventory) to ensure adequate respirator availability and product usage within these optimized shelf lives.

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LIST OF ABBREVIATIONS

CADTH	Canadian Agency for Drugs and Technologies in Health
CDC	Centres for Disease Control and Prevention, US
FDA	Food and Drug, US
FFR	Filtering facepiece respirator
INESSS	Institut national d'excellence en santé et services sociaux, QC
PPE	Personal protective equipment
NICE	National Institute for Health and Care Excellence, UK
NIOSH	National Institute for Occupational Safety and Health, US
WHO	World Health Organization

USE OF N95 RESPIRATORS PAST THEIR SHELF LIFE

1. BACKGROUND

The COVID-19 pandemic necessitated the stockpiling of personal protective equipment (PPE) as a way to ensure the availability of necessary quantities during the public health emergency. Although The National Institute for Occupational Safety and Health (NIOSH) does not require manufacturers of filtering facepiece respirators (FFR) to assign an expiry date for certification purposes, most FFR manufacturers designate a 5-year shelf life. During and following the pandemic, healthcare establishments developed stockpiled reserves of PPE, including N95 respirators, that have since exceeded their manufacture-designated shelf life.

The main objective of N95 masks is to protect healthcare personnel from infection with airborne pathogens by preventing aerosol transmission, which raises important questions about the continued effectiveness of expired masks. Given the considerable procurement and storage costs of PPE and the potential need for such equipment in future public health crises, the destruction of large quantities of potentially effective respirators may not be fiscally responsible.

It is therefore important to evaluate whether N95 respirators past their shelf life continue to be safe and effective in protecting healthcare personnel.

1.1 Reason for Brief Report

The MUHC currently has a stockpile of 286,100 N95 Champak F550 respirators that have recently expired or are due to expire in 2024. These masks were purchased by the MUHC at a cost of \$228,654.

This stockpile resulted due to a confluence of factors:

- A transition from Champak F550 respirators to Medicom masks in 2022
- Offloading of Champak F550 respirators from the MSSS to deplete their inventory

The Champak F550 masks are NIOSH approved, were deemed safe by an evaluation of the Occupation Health and Safety department of the MUHC, and underwent extensive fit testing. Given the considerable acquisition, storage and disposal costs, as well as potential safety risks to healthcare personnel, there is a need to evaluate the risks and benefits of using N95 respirators past their expiration date.

Therefore, this evaluation was requested by Kim Tanguay, Nursing Practice Consultant, Directorate of Nursing at the McGill University Health Centre (MUHC) on November 10, 2023.

2. POLICY AND EVALUATION QUESTIONS

2.1 Policy Question

Can the MUHC continue to use N95 respirators past their expiration date?

2.2 Evaluation Questions (Objective of this report)

- What is the evidence on the safety and effectiveness of expired N95 respirators in terms of aerosol filtration efficiencies?
- Is there guidance for the safe use of N95 respirators past their expiry date?

3. METHODS

3.1 Literature search

We conducted a literature search of PubMed, health technology assessment (HTA) databases (CADTH, INESSS, NICE) and regulatory bodies (WHO, CDC, Health Canada) for studies on effectiveness of and guidance for expired N95 respirators. The most recent search was conducted on November 17, 2023. The following key words were used in PubMed: (N95) AND ((expired) OR (stockpile) OR (shelf life)).

4. EVIDENCE FROM SCIENTIFIC LITERATURE

4.1 Results of the literature search

We identified three studies that evaluated the filtration performance of expired N95 respirators, one framework for evaluating stockpiled N95 respirators, and several guidance documents from regulatory bodies relating to the use of expired respirators in times of low supply.

4.2 Factors that determine N95 respirator performance

4.2.1 Aerosol filtration efficiencies

Filtration performance of N95 respirators is based on NIOSH Standard Test Procedures. The Pass/Fail criteria for Particulate Filter Efficiency for N95 is $\geq 95\%$, i.e. **N95 FFRs must have a particle penetration of less than 5.0%.**

- The largest study to date to evaluate the filtration performance of stockpiled N95 respirators was conducted by **NIOSH** between 2017 and 2019 to understand the viability of respirators that have undergone long-term storage.¹ They evaluated 11 types of N95 filtering facepiece respirators manufactured between 2003 and 2013, for a total of 3,971 respirators sampled from 91 manufacturing lots. The models included were from the following manufacturers: 3M, Gerson, Kimberly Clark, Medline/Alpha Protech and Moldex. The expired respirators were compared to a control group of new respirators of the same model and evaluated for inhalation and exhalation resistance and filtration performance, based on NIOSH Standard Test Procedures.

Findings: NIOSH found that 98% (3,895 out of the 3,971) of the evaluated respirators maintained their inhalation and exhalation resistance and filtration performance, with the upper limit of the 99% confidence intervals well within the NIOSH minimum performance requirements.

- In a small quality improvement study by **Sickbert-Bennett et al.** in 2020 that tested 29 fitted face mask models, the authors report that N95 respirators from the company 3M that were 9 and 11 years past their expiry date had fitted filtration efficiencies of 98.5% and 97%, respectively.²

Findings: The authors conclude that expired N95 respirators with intact elastic bands continue to provide effective protection against submicron particle penetration. However, this conclusion was drawn on the basis of testing only 3 masks of each model (3M 8210 and 3M 1816) on a single male participant.

- An older study by **Viscusi et al.** in 2009 evaluated the filtration performance of 63 masks from 21 different N95 respirator models that had been stored in their original packaging for an average of 7.3 years at temperatures ranging between 15°C and 32°C and relative humidity between 20% and 80%.³

Findings: 19 of 21 models passed the NIOSH Standard Test Procedure and had particle penetration of less than 5.0%; 2 models failed the test by less than 1.0%.

The authors report that there did not appear to be any correlation between the length of storage and failure to pass the initial penetration test, with models stored for approximately 10 years still achieving very good filtration performance.

4.2.2 Respirator strap performance

A study by **Rottach et al.** examined the degradation of polyisoprene (synthetic rubber) elastomer straps taken from a sample of two N95 FFR models (n=54) of various ages (1-10 years).⁴ The masks were stored in poor conditions, and factors like high heat can deform strap stability.

Findings: The authors report that the tension of the straps from one model were found to differ according to age, while the second model showed no correlation with time. They also note that minor variations in strap tensile properties may not result in significant differences in respirator seal quality. Based on this small study in 2 respirator models, they conclude that prolonged storage may affect the tensile properties of headstraps for some models of N95.

4.2.3 Storage conditions

There is consensus from regulatory bodies that the following storage conditions could affect the protective ability of N95 respirators:⁵

- Length of storage time beyond its shelf life;
- Storage conditions such as temperature, humidity and exposure to dust, moisture, sunlight and direct light, proximity to fans, windows, doors, and ventilation systems, and presence of mould;¹
- Damage to pallets, cases, boxes and individual masks, which could result in damage of nose bridge, nose foam, filter media or straps;¹
- Tensile properties of straps after prolonged storage.

4.3 Summary of NIOSH findings

- NIOSH acknowledges that shelf life issues will be an ongoing consideration for managing stockpiled respirators.¹
- Results of their large study provide important information on the levels of protection that may be provided by respirators stored for long periods in U.S.

stockpiles. Nonetheless, the results may not be applicable to other respirator models or stockpile facilities.¹

- They report that the majority of respirators evaluated would be protective across a variety of storage conditions, but important variability was observed.
- NIOSH does not have enough information to definitively know the level of protection that may be provided by all respirators that 1) are stored for prolonged periods of times; 2) are stored under various storage conditions; or 3) have exceeded the approval holder's designated shelf life.¹

5. SOLUTIONS FOR STOCKPILED RESPIRATORS

5.1 NIOSH recommendations for stockpiled respirators

- Given the current reality that most US facilities have stockpiles of respirators and other PPE stored past their designated shelf life, there is need for uniform guidance related to inventory management for expired respirators.¹
- **One solution is the development of an economically sustainable stockpile quality assurance process, such that stockpiles from across the nation test their respirator inventories using NIOSH's minimum performance requirements, and share their results with product manufacturers to optimize shelf life designations.**¹
- In addition to optimizing shelf life designations, NIOSH recommends a national stockpile maintenance approach to purchase, store, and deploy respirators, for example by rotating inventory, to ensure adequate respirator availability and usage within these optimized shelf lives.¹

5.2 A framework for evaluating and testing stockpiled FFRs

Yoon et al. developed a 4-step process aimed at stockpile managers and policy-makers for evaluating the performance viability of N95 FFRs:⁶

1. Define the population of FFRs to be sampled: "If quality information is desired for a single FFR model with a few lots, designing a study that allows for quality inferences to be made to the entire FFR model (i.e., the study population) could be easily accomplished by randomly sampling FFRs from each lot."

2. Sampling strategies: Once the population is defined, a random sampling process should be used to ensure that each unit within a lot has an equal probability of being selected for testing. When lot information is available and large quantities of FFRs are stored, testing and evaluation may use the Lot Quality Assurance Sampling (LQAS) technique. This is a stratified random sampling design developed to facilitate decisions to accept or reject the entire lot, by allowing generalization of the results from sample data to the wider population from which the sample was drawn (Appendix [Figure 1](#))
3. Inspecting and testing the sample: Samples should be tested for:
 - Visual inspection: For visible damage, deformation or presence of dust or mould (see example visual checklist in Appendix [Figure 2](#))
 - Standard Testing Procedures (STP): STPs should be evaluated in accordance with government or Occupational Health and Safety norms and regulations. NIOSH STP-0059 to determine Particulate Filter Efficiency Level should be considered the primary performance test to evaluate the viability of the FFRs being tested. For N95 FFRs, particulate penetration is required to be at or below 5.0%.
4. Evaluating the test results: The number of units per sample that fail the above criteria can be tabulated and then extrapolated to the entire population to determine the probability of observing a similar failure rate in the whole population of FFRs.

6. GUIDANCE FOR EXPIRED N95 RESPIRATORS

Most regulatory bodies including Health Canada and the WHO agree that respirators beyond their labelled shelf life, while no longer considered approved, may retain adequate filter performance if stored properly.⁵ However, because loss of elasticity in the bands or straps can impact a proper seal, they recommend that respirators be closely inspected for signs of damage and fit-tested before use.

6.1 Health Canada guidance⁷

- An expired mask can still be effective at protecting the health care provider if:
 - The straps are intact
 - There are no visible signs of damage

- They can be fit-tested
- Health care providers should inspect the mask and perform a seal check
- There is no specific timeframe beyond the expiry dates for N95 respirators at which they would no longer be considered suitable for use.

6.2 WHO guidance

The WHO issued guidance on the use of PPE beyond the manufacturer-designated shelf life or expiration date in times of PPE shortages.⁸ Although respirators past their designated shelf life are no longer considered approved in accordance with their associated regional/international standards, the WHO [recommends](#) that:

- Such items may still be effective and should be inspected before use to verify equipment is in good condition with no degradation that could affect performance. Specifically, a respirator may still be effective if
 - It has been appropriately stored to avoid the effects of moisture or contamination,
 - The straps have remained intact,
 - There are no visible signs of damage and
 - A self-fit test/seal check can be performed successfully by the wearer before use.

7. CONCLUSIONS

- The best available evidence on the performance of expired or stockpiled N95 respirators is derived from a NIOSH evaluation of 3,971 stockpiled N95 respirators manufactured between 2003 and 2013, which found that 98% of the evaluated respirators maintained their inhalation and exhalation resistance and filtration performance after their expiry date.
- Nonetheless, NIOSH cautions these results may not be applicable to other respirator models not included in this study or stockpile facilities with various storage conditions. The Chompak F550 N95 respirators were not included in this study.
- Most regulatory bodies including Health Canada, the FDA and the WHO agree that respirators beyond their labelled shelf life, while no longer considered approved,

may retain adequate filter performance if stored properly. However, because loss of elasticity in the bands or straps can impact a proper seal, they recommend that respirators be closely inspected for signs of damage and fit-tested before use.

- To guide stockpile personnel, NIOSH recommends the creation of a quality assurance process for stockpiled respirators to determine optimal shelf life designations.
- NIOSH has also developed a 4-step framework to guide stockpile managers in evaluating the performance of stockpiled respirators, such that a representative sample of respirators can undergo quality control to determine the probability of observing a similar failure rate in the whole population of FFRs.

8. KEY FINDINGS AND RECOMMENDATIONS

- There is reasonable evidence to indicate that N95 respirators past their designated shelf life retain adequate filter performance if stored properly. NIOSH reported that 98% of almost 4000 evaluated N95 respirators past their expiry date maintained their protective properties. However, these results may not be generalizable to respirators not included in their analysis, including Champak F550 N95 respirators.

Recommendations for the MUHC:

- Stockpile managers should perform quality assurance testing on representative samples of Champak F550 N95 respirators from different lots and storage conditions, as described by NIOSH, so the viability of the complete stockpile can be ascertained. Representative sample testing allows extrapolation of findings to the population of respirators from which the sample was derived, allowing for a cost-effective solution;
- The following factors should be evaluated:
 - Visual signs of damage to respirators and fit-testing before use
 - Filtration efficiency
 - Storage conditions

Recommendations at the national/province level (NIOSH recommendation):

- The creation of **an economically sustainable stockpile quality assurance process**, to ensure that stockpiles of respirator inventories across the province are

evaluated using NIOSH's minimum performance requirements. Results should be shared with product manufacturers to optimize shelf life designations.

- A stockpile maintenance approach at the national level to purchase, store, and deploy respirators (e.g., by rotating inventory) to ensure adequate respirator availability and product usage within these optimized shelf lives.

REFERENCES

1. Greenawald L, Moore S, Yorio PL. Inhalation and Exhalation Resistance and Filtration Performance of Stockpiled Air-Purifying Respirators: Overall Performance of Nearly 4,000 Respirators Sampled from Ten Stockpile Facilities. Pittsburg, PA U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, National Institute for Occupational Safety and Health. 2020. Available from: <https://www.cdc.gov/niosh/npptl/ppecase/pdfs/PPE-CASE-Aggregated-Stockpile-Study-03252020-508.pdf>.
2. Sickbert-Bennett EE, Samet JM, Clapp PW, et al. Filtration Efficiency of Hospital Face Mask Alternatives Available for Use During the COVID-19 Pandemic. *JAMA Intern Med.* Dec 1 2020;180(12):1607-1612.
3. Viscusi DJ, Bergman M, Sinkule E, Shaffer RE. Evaluation of the filtration performance of 21 N95 filtering face piece respirators after prolonged storage. *Am J Infect Control.* Jun 2009;37(5):381-386.
4. Rottach DR, Lei Z. Stockpiled N95 Filtering Facepiece Respirator Polyisoprene Strap Performance. *J Int Soc Respir Prot.* 2017;34(2):69-80.
5. Health Quality Ontario. Extended Use and Layering of N95 Respirators and Use of Expired Personal Protective Equipment. Canada. 2020 Available from: <https://www.hqontario.ca/Portals/0/documents/evidence/reports/extended-use-and-layering-of-n95-respirators-and-use-of-expired-ppe.pdf>.
6. Yoon KN, Greenawald LA, Rottach DR, Pollard JP, Yorio PL. A General Framework to Test and Evaluate Filtering Facepiece Respirators Considered for Crisis Capacity Use as a Strategy to Optimize Supply. *J Int Soc Respir Prot.* 2020;36(1):36-51.
7. Health Canada. Optimizing the use of masks and respirators during the COVID-19 outbreak [Internet]. . Ottawa, ON: Health Canada; 2020. Available from: <https://www.canada.ca/en/health-canada/services/drugs-health-products/covid19-industry/medical-devices/personal-protective-equipment/medical-masks-respirators/health-professionals.html - a32>.
8. World Health Organization. Rational use of personal protective equipment for COVID-19 and considerations during severe shortages. 2020. Available from; https://iris.who.int/bitstream/handle/10665/338033/WHO-2019-nCoV-IPC_PPE_use-2020.4-eng.pdf?isAllowed=y&sequence=1.

APPENDIX


Figure 1: Hypothetical stratified random sampling design for N95 FFRs

Source: Yoon, K. N., Greenawald, L. A., Rottach, D. R., Pollard, J. P., & Yorio, P. L. (2020). A General Framework to Test and Evaluate Filtering Facepiece Respirators Considered for Crisis Capacity Use as a Strategy to Optimize Supply. *Journal of the International Society for Respiratory Protection*, 36(1), 36–51.

Visual Inspection Items	Evaluation
1. Does the packaging have any signs of damage ?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA
2. Does the respirator have an odor ?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA
3. Does the respirator have a presence of dust ?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA
4. Is the respirator deformed in any way ?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA
5. Is the nose clip cracked, corroded, or detached?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA
6. Is the nose foam flaking, or does it appear damaged?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA
7. Does the respirator appear moldy ?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA
8. If applicable, are the staples attaching the straps cracked, corroded, or rusty ?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA
9. Are the straps detached from the respirator?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA
10. Are the straps damaged in any way?	<input type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> NA

Figure 2. Example Respirator Visual Inspection Checklist

Source: Yoon, K. N., Greenawald, L. A., Rottach, D. R., Pollard, J. P., & Yorio, P. L. (2020). A General Framework to Test and Evaluate Filtering Facepiece Respirators Considered for Crisis Capacity Use as a Strategy to Optimize Supply. *Journal of the International Society for Respiratory Protection*, 36(1), 36–51.