NOTE DE SERVICE
MEMORANDUM

OPTILAB Montréal-CUSM

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==========

Dear Colleagues

Beginning **November 3, 2026**, and gradually, the Microbiology Laboratory at the Glen site will change the way respiratory viruses are detected in clinical samples. The laboratory will now use the **Cobas Respiratory Flex assay**, a newly available multiplex assay approved by Health Canada.

This assay detects **15 different respiratory virus species**, including **SARS-CoV-2**, eliminating the need for supplementary testing for SARS-CoV-2.

The additional viral targets include:

- Influenza A
- Influenza B

- Respiratory Syncytial Virus (RSV)
- Human Metapneumovirus (HMPV)
- Adenovirus
- Enterovirus/Rhinovirus
- Parainfluenza Viruses 1, 2, 3, and 4
- Human Coronaviruses HCoV-OC43, HCoV-229E, HCoV-NL63, and HCoV-HKU1 (reported collectively as **"common coronaviruses"**)

The new assay will replace several currently available tests, including:

In-house multiplex respiratory virus assay	Respiratory BioFire Panel	Cobas SARS-CoV-2
Cobas SARS-CoV-2 / Influenza A / Influenza B	GeneXpert SARS-CoV-2	GeneXpert Quad (SRAS-CoV-2, Influenza A, Influenza B et VRS)
Simplexa SARS-CoV-2		

These assays will no longer be available. Requests submitted for any of these tests will be converted by the laboratory to the new assay.

Ordering the Assay

The test can be ordered in **OACIS** as:

Respiratory Viruses & SARS-CoV-2
(*"MultiResp & SCoV-2"*)

Required Specimens: Nasopharyngeal (NP) Swabs

- NP swabs collected in **Universal Transport Medium (UTM)**, such as **Copan UTM**.

Important: NP swabs submitted in **COBAS PCR medium** (yellow-cap tubes) cannot be processed with the current version of the assay.During a one-month transition period, NP swabs received in COBAS medium will continue to be analyzed using existing methods (the in-house multiplex assay and Cobas SARS-CoV-2 assay). After the transition period, specimens received in COBAS medium will be cancelled and a request for a new specimen will be issued.

Other Respiratory Specimens

The following specimen types have not yet been fully validated:

- Bronchoalveolar lavage (BAL) specimens
- Endotracheal aspirates
- Nasopharyngeal aspirates

These specimens will be tested; however, results will be reported with a cautionary note regarding interpretation.

Turnaround Time

Turnaround time will vary depending on the specimen's arrival time in the laboratory relative to the next scheduled testing run. However, in all cases, results are expected to be available within **12 hours**.